

Laminar Flow Hood

Laminar cleanliness, safeguarding critical processes



LAMINAR FLOW HOOD

Laminar cleanliness, safeguarding critical processes


A laminar flow hood is an air purification device designed to provide a highly clean environment for localized areas. It primarily consists of an enclosure, fan, high efficiency particulate air (HEPA) filter, diffuser plate, lighting systems, and other components. The unit creates a stable unidirectional airflow (vertical/horizontal) by forcing air through HEPA or ULPA filters, effectively isolating the work zone from external contaminants while maintaining superior cleanliness standards.

This product features a compact structure and easy operation, and can be optionally equipped with lighting systems, differential pressure gauges, pre-filters, and other functional components. It is suitable for critical processes with stringent requirements, such as aseptic operations, sample protection, and precision assembly.



- 1 Air inlet grille
- 2 DOP/PAO test port
- 3 Differential pressure gauge
- 4 LCD display
- 5 Airflow equalizing membrane/diffuser plate

Working Principle



1. Air Intake

The laminar flow hood draws in ambient air through the top or rear air inlet grille. In some models, the air first passes through a coarse air filter to remove large dust particles.

2. Filtration & Purification



The incoming air passes through a high-efficiency particulate air (HEPA) filter (efficiency $\geq 99.99\%$ @ $0.3 \mu\text{m}$) or an ultra-low penetration air (ULPA) filter (efficiency $\geq 99.9995\%$ @ $0.12 \mu\text{m}$), filtering out the vast majority of particulates and microorganisms.

3. Formation of Laminar (Unidirectional) Flow

Filtered air is delivered at a constant velocity and in a constant direction from beneath the filter (vertical laminar flow) or from the front (horizontal laminar flow), forming a stable laminar airflow that pushes contaminants outward from the work zone.

4. Contaminant Removal

Particles or contaminants generated during operations are quickly swept away by the airflow instead of retaining in the work zone, thereby maintaining area cleanliness.



Specifications



Material

The enclosure and all parts are constructed from SUS304.



Pressure

Negative or positive pressure.



Average Air Velocity

0.45±20% m/s.



Illuminance

≥ 300 LUX.



Functions

Air velocity feedback, variable frequency speed control, multi-parameter alarm system, and touchscreen operation.



Air Inlet

Side or top entry.



Cleanliness Class

Class A.



Noise

≤ 65 dB.



Testing

PAO scan leak test port.



Model	Rated Airflow	External Dimensions (L × W × H) mm	Work Zone Dimensions mm	HEPA Specification	Power Supply
	m ³ /h			mm	
JLFZJ-600	600	730 × 730 × 790	1220 × 610	610 × 610 × 90 1pcs	AC 220V/50Hz
JLFZJ-1200	1200	1350 × 730 × 790	1220 × 610	610 × 610 × 90 2pcs	
JLFZJ-1800	1800	1350 × 1030 × 790	1220 × 915	915 × 610 × 90 2pcs	
JLFZJ-2400	2400	1350 × 1350 × 790	1220 × 1220	1220 × 610 × 90 2pcs	
JLFZJ-4800	4800	2830 × 1030 × 790	2440 × 915	915 × 610 × 90 4pcs	

Application

Widely used in industries requiring stringent air cleanliness standards, including pharmaceuticals, biological laboratories, electronics manufacturing, and precision instrument assembly.



PureFlow

Hengshui Pureflow Tech Co., Ltd.

- +86-13091189390 (WhatsApp & WeChat)
- hugo@pureflowtec.com
- www.pureflowtec.com



WeChat



WhatsApp